

LightMix[®] Kit *Mycoplasma hominis* (EC)

Cat.-No. 40-0139-32

Internal Control (IC) changed to Extraction Control (EC), 32 rxn/vial

Kit with reagents for the detection of *Mycoplasma hominis* DNA using the Roche Diagnostics LightCycler[®] 1.x / 2.0 / 480 II Instruments or cobas z 480 Analyzer.

Lyophilized mix of primers and probes for a total of 96 reactions with a final volume of 20 µl each.

Store protected from light at room temperature (18-25°C), do NOT freeze!

Instructions for use with the LightCycler[®] 1.x / 2.0 Instruments see pages 4-5

Instructions for use with the LightCycler[®] 480 II and cobas z 480 Analyzer see pages 6-7

1. Introduction

Mycoplasma are the smallest bacteria known. Since they lack a cell wall they are unassailable for the common beta-lactam antibiotics targeting cell wall synthesis, however, they are sensible to Tetracyclines. *Mycoplasma hominis* (*M. hominis*) is commensal of the gastrointestinal tract and not harmful for immunocompetent individuals. Infections of the urogenital tract in adults are frequent and mostly symptomless, but *M. hominis* is suspected to cause pelvic inflammatory disease and even infertility. In premature infants *M. hominis* can cause sepsis or bacterial meningitis resistant to antibiotic treatment.

The glyceraldehyde-3-phosphate dehydrogenase (*gap*) gene is a common target for PCR based tests¹.

2. Description

This kit provides a fast and accurate system to detect *M. hominis* genome DNA in a nucleic acid extract. A 129 bp long fragment of the glyceraldehyde-3-phosphate dehydrogenase (*gap*) gene from *M. hominis* is amplified with specific primers and detected with probes labeled with LightCycler[®] Red 640.

The control reaction generates a 125 bp fragment from the PhHV target, detected with LightCycler[®] Red 690 labeled hybridization probes. This PCR has no visible impact on the *M. hominis* specific reaction and will even fail in the presence of higher amounts of target (1,000 copies and more).

The former internal control (IC) has been changed to a spiked extraction control (sEC) to monitor a successful extraction and demonstrate the ability to run an amplification reaction (no PCR inhibition). We recommend to use the 'Extraction Control' procedure; for compatibility with the former procedure the use as IC is also described. The extraction control target ¹ECT contains Lambda and PhHV DNA and may be used also for other LightMix kits, without adding the ECT provided with the other kit(s).

Note: With a MagNA Pure Compact extraction the recovery rate for the EC target can be very low or the DNA target even gets lost; the amount of ECT might have to be adapted to the extraction method.

This kit can be combined with the LightMix[®] Kit *M. genitalium* (40-0169-32). The species present in the samples can be identified through the specific melting points at 67-69°C for *M. genitalium* and 56-64°C for *M. hominis*. A laboratory protocol for combination of both assays is described in this manual under **6.3 Preparation of the LightCycler[®] reaction mix**. The control reaction used for both kits is identical.

Laboratories running LightCycler[®] 2.0 / 480 II or cobas z 480 instruments may switch to use the 'LightMix[®] Kit 40-0460-32 M. hom/gen U.urea EC' to detect the pathogens *M. genitalium*, *M. hominis* and *U. urealyticum* combined with a spiked Extraction Control in one single multiplex reaction.

The use of a color compensation file generated with the ColorCompensation kit 40-0318 is a prerequisite to run a prerequisite to detect the control reaction.

The supplied standard row allows to determine the linear range of the reaction and to estimate the quantity of the target sequence in unknown samples.

The kit must be used with 'LightCycler[®] FastStart DNA Master HybProbe' only (capillary and plate based LightCycler[®] Instruments).

3. Set Contents

- 3 Vials with **green** cap containing premixed primers and probes for each 32 reactions of *M. hom.*
- 3 Vials with **white** cap containing premixed primers and probes for each 32 control reactions
- 1 Standard row with 6 lyophilized standards *M. hominis* 10¹ to 10⁶ target equivalents per rxn
- 1 Sealing foil for the standard row
- 1 Vial with **white** cap containing Extraction Control Target (ⁿECT): 4.8 x 10⁶ copies (total)
- 1 Vial with **black** cap containing the stabilizer for the No Template Control (NTC)
- 1 Certificate of Analysis (CoA) with lot-specific data

4. Additional Reagents and Items Required

LightMix [®] Kit ColorCompensation HybProbe 40-0318-00	Roche Diagnostics Cat.-No. 05 997 704 001
LightCycler [®] FastStart DNA Master HybProbe	Cat.-No. 03 003 248 001
LightCycler [®] Capillaries (20 µl) (LightCycler [®] 1.x / 2.0 Instruments)	Cat.-No. 04 929 292 001
LightCycler [®] 480 Multiwell Plate 384, white (LightCycler [®] 480 Instrument)	Cat.-No. 04 729 749 001
or LightCycler [®] 480 Multiwell Plate 96, white (LightCycler [®] 480 Instrument)	Cat.-No. 04 729 692 001

For use in LightCycler[®] 1.x Instruments with software version 3.5.3 read channel F2 instead of 640 and F3 instead of channel 705. We recommend to upgrade to software version 4.10 or higher.

4.1. Optional Additional Reagents

High Pure PCR Template Preparation Kit	Cat.-No. 11 796 828 001
--	-------------------------

5. Product Characteristics

PCR results (activation, 50 cycles and melting curve) are obtained within 50 minutes with the capillary based LightCycler[®] 1.x / 2.0 Instruments and within 80 minutes with '480' plate based Instruments.

Sensitivity

Reagents detect 10 copies/reaction of positive control target DNA using FastStart DNA Master HybProbe.

Measuring range

The linear measuring range of the assay is 10² to 10⁶ copies of *Mycoplasma hominis* DNA.

Storage and Stability

- Lyophilized reagents are stable for at least 6 months after shipment when stored protected from light at room temperature (18-25°C). See expiry date on the product label.
- **Do not freeze** lyophilized reagents.
- Dissolved reagents are stable for at least 10 days when stored protected from light and refrigerated (4°C).

6. Experimental Protocol

Start programming before preparing the solutions. See the Instrument operator's manual for details.

6.1. Preparation of Parameter-Specific Reagents and reagents for the EC (32 reactions):

One reagent vial with a **green** cap contains primers and probes to run 32 reactions of *M. hominis*.
One reagent vial with a **white** cap contains primers and probes to run 32 control reactions.

Check for the colored pellet, then **add 66 µl** PCR-grade water, mix (vortex) and spin down.

► Use **2 µl** reagent for a 20 µl PCR reaction.

| This solution is stable at least ten days when stored refrigerated at 4°C. Avoid prolonged exposure to light.

6.2. Preparation of Extraction Control Target (ECT):

Add **1,200 µl** PCR-grade water to the vial with **white** cap containing the Extraction Control Target

6.3. Sample Preparation:

Before extraction add **10 µl** of ECT (vial **white** cap) to the sample; the amount may have to be adapted to the extraction method to get a Cp value in the range of 28-32. **Skip if IC procedure is used.**
Perform nucleic acid preparation as described in the protocol of the extraction kit used (see 4.1).

6.4. Preparation of No Template Extraction (NTC):

Add **900 µl** PCR-grade water to the **NTC** (vial **black** cap) and add **100 µl** of ECT (vial **white** cap).

► Use **5 µl** No Template Control DNA for a 20 µl PCR reaction.

6.5. Preparation of the Standard Row

The target DNA is provided in 6 different quantities to yield from 10 to 10⁶ target molecules in 5 µl once resolved. Start with the lowest concentrated standard (first tube from the extended lip). Use the pipette tip to punch a hole through the sealing foil. Add **40 µl** of **NTC** (vial **black** cap) to each well. **Use water if the IC procedure is selected.** Mix the target DNA by pipetting the solution up and down 10 times.



► Use **5 µl Standard** for a 20 µl PCR reaction

| This standard solution is not long-term stable. Use only fresh prepared solutions only. After adding the target DNA to the reaction mix, use the provided sealing foil to close the vials in order to avoid contaminations.

6.6. Preparation of the Reaction Mix

Include Positive Controls and at least one 'No Template Control' (NTC). In a cooled tube, prepare the reaction mix by multiplying the single reaction volumes by the number of reactions plus one reserve :

sEC procedure		For use with the Roche FastStart Master	Optional: IC	
Single reaction	Duplex reaction	Component	Single reaction	Duplex reaction
6.6 µl	4.6 µl	Water, PCR-grade (colorless cap, provided with the Roche Master kit)	6.1 µl	4.1 µl
2.4 µl	2.4 µl	Mg ²⁺ solution 25 mM (blue cap, provided with the Roche FastStart kit)	2.4 µl	2.4 µl
2.0 µl	2.0 µl	PSR mix (parameter specific reagents, see 6.1 for <i>M. hominis</i>)	2.0 µl	2.0 µl
-	2.0 µl	PSR mix (parameter specific reagents, see 6.1 for <i>M. genitalium</i>)	-	2.0 µl
2.0 µl	2.0 µl	Primers and probe mix for the IC/EC	2.0 µl	2.0 µl
-	-	ECT (white cap, DNA control target, see 6.2.)	0.5 µl	0.5 µl
2.0 µl	2.0 µl	Roche Master (red cap, for preparation see Roche manual)	2.0 µl	2.0 µl

15.0 µl

Volume of reaction mix

15.0 µl

Table 1

To run the assay without the control reaction substitute ECT with 0.5 µl PCR-grade water.

Mix gently, spin down and **transfer 15 µl** of the reaction mix to a capillary or well.

Add 5 µl of sample or standard to each capillary or well for a final reaction volume of 20 µl.
Close the capillaries / attach a foil to the multiwell plate and seal, and spin down.

Start run.

7. LightCycler® 1.x / 2.0 Instruments

7.1. Programming

The protocol consists of four program steps

- 1: Denaturation: sample denaturation and enzyme activation
- 2: Cycling: PCR-amplification of the target DNA
- 3: Melting: melting curve analysis for identification of the PCR product derived from the target DNA
- 4: Cooling: cooling the instrument

Program Step:	Denaturation	Cycling			Melting			Cooling
Parameter								
Analysis Mode	None	Quantification mode			Melting Curves mode			None
Cycles	1	50			1			1
Target [°C]	95	95	62	72	95	40	85	40
Hold [hh:mm:ss]	00:10:00	00:00:05	00:00:05	00:00:15	00:00:20	00:00:20	00:00:00	00:00:30
Ramp Rate [°C/s]	20	20	20	20	20	20	0.2	20
Sec Target [°C]	-	-	55	-	-	-	-	-
Step Size [°C]	-	-	0.5	-	-	-	-	-
Step Delay (Cycles)	-	-	1	-	-	-	-	-
Acquisition Mode	None	None	Single	None	None	None	Cont.	None

Table 2

7.2. Data Analysis

Switch the color compensation mode on. If this mode is not enabled run the color compensation program. Follow the instructions in the manual of 'LightMix® Kit – Color Compensation HybProbe.

Perform data analysis, as described in the LightCycler® Instrument operator's manual.

We recommend using the Second Derivative Maximum method (Automated (F" max)). The cycle number of the Crossing Point (Cp) of each sample is calculated automatically. The Fit Points method is more-error prone due to the user's influence.

View *M. hominis* data in channel 640, Quantification mode. The negative control (NTC) must show no signal.

For the identification of the PCR product view *M. hominis* (and *M. genitalium* when combined) data in channel 640, Melting Curves mode: specific melting points at 60-64°C (eventually 56°C) for *M. hominis* (and 67-69°C for *M. genitalium*, if assays are combined in one reaction).

For the Control Reaction view channel 705 data. The negative control and the low-concentrated *M. hominis* DNA samples (10 to 1,000 copies) should show an amplification curve for the EC/IC with a Cp at approximately cycles 30-32.

The provided standard row of cloned DNA with concentrations in the range from 10⁶ copies/rxn to 10¹ copies/rxn of *M. hominis* should have Cp values between cycles 18 and 37.

7.3. Interpretation of Data

Sample 640 <i>M. hominis</i>	Sample 705 Ctrl Reaction	Channel 640 PositiveControl	Negative Control (NTC)	Result (warnings)
no amplification	detectable	amplification	negative	Negative (not detectable)
Cp < 39 ⁺	not relevant	amplification	negative	Positive for <i>M. hominis</i>
no amplification	not detectable	amplification	not relevant	PCR failure, repeat experiment
not relevant	not relevant	no amplification	not relevant	PCR failure, repeat experiment
not relevant	not relevant	not relevant	positive	Contamination, repeat experiment

Table 3. Typical analysis results (LightCycler® 2.0 Instrument, Roche Master: FastStart)

⁺ The cut off value is a recommendation only - this value has to be defined by the user. Compare with the lot specific values for 10 copies as reported in the Certificate of Analysis (CoA). 1-2 cycles later than observed for 10 copies corresponds to ~5 copies.

7.4. Sample Data – Typical Results




Fig.1. LightCycler® 2.0 sample data for the *M. hominis* detection system.

Upper panels: Left panel channel 640 quantification mode (Sec. Derivative Maximum) with amplification curves for *M. hominis*. Right panel channel 640 melting analysis for *M. hominis*.

Central panels: Left panel channel 705 quantification mode (Second Derivative Maximum) for the control reaction. Right panel channel 705 melting analysis for the control reaction (not relevant for detection).

Lower panels: Left panel channel 640 quantification mode (Second Derivative Maximum) for *M. hominis* and *M. genitalium*. Right panel channel 640 melting analysis for *M. hominis* and *M. genitalium* - identification of the species if combined with LightMix® Kit 40-0169 *M. genitalium*.

Note: Fluorescence levels may vary also depending on instrument settings. The characteristics of the curve (shape and trend of fluorescence levels) should be similar to the curve in the manual. The fluorescent curves over cycles (quantification mode) must be smooth and not zigzag shaped. The Cp values will vary from instrument to instrument by up to 2 cycles, while the distance between two dilution steps should be relatively constant (delta Cp). The Cp values described in this manual (chart text) have been obtained using the supplied standard row; the entire set of curves can be horizontally shifted.

8. LightCycler® 480 II Instrument and cobas z 480 Analyzer

8.1. Programming

The protocol consists of four program steps:

- 1: Denaturation: sample denaturation and enzyme activation
- 2: Cycling: PCR-amplification of the target DNA
- 3: Melting: melting curve analysis for identification of the PCR product derived from the target DNA
- 4: Cooling: cooling the instrument

Detection Format:

LightCycler® 480 II Instrument: 465-510, 498-640, 498-660

cobas z 480 Analyzer: 465-510, 498-645, 498-700

Program Step:	Denaturation	Cycling			Melting			Cooling
Parameter								
Analysis Mode	None	Quantification mode			Melting Curves mode			None
Cycles	1	50			1			1
Target [°C]	95	95	62	72	95	40	85	40
Hold [hh:mm:ss]	00:10:00	00:00:05	00:00:05	00:00:15	00:00:30	00:02:00	00:00:00	00:00:30
Ramp Rate [°C/s] 96	4.4	4.4	2.2	4.4	4.4	1.5	-	1.5
Ramp Rate [°C/s] 384	4.6	4.6	2.4	4.6	4.6	2.0	-	2.0
Sec Target [°C]	-	-	55	-	-	-	-	-
Step Size [°C]	-	-	0.5	-	-	-	-	-
Step Delay (Cycles)	-	-	1	-	-	-	-	-
Acquisition Mode	None	None	Single	None	None	None	Continuou	None
Acquisitions [per °C]	-	-	-	-	-	-	1	-

Table 4

8.2. Data Analysis

Note: cobas z 480 Analyzer signal levels are about 50% compared to LightCycler® 480 II results.

Switch the color compensation mode on. If this mode is not enabled run the color compensation program. Follow the instructions in the manual of the 'LightMix® Kit Color Compensation HybProbe.

Perform data analysis, as described in the Instrument operator's manual. We recommend using the Second Derivative Maximum method (Automated (F" max). The cycle number of the Crossing Point (Cp) of each sample is calculated automatically.

View data with Filter Combination 498-640, Quantification mode. The negative control (NTC) must show no signal. For identification of the PCR product view data with Filter Comb. 498-640, Melting Curves mode with a Tm of 56-64°C for *M. hominis* (and 67-69°C for *M. genitalium*, if combined).

If the control reaction is used, view data with Filter Combination 498-660, Quantification mode. The negative control and the low-concentrated *M. hominis* DNA samples (10 to 1,000 copies) should show an amplification curve for the EC/IC with a Cp at approximately cycle 28-30.

The provided standard row of cloned DNA with concentrations in the range from 10⁶ copies/rxn to 10 copies/rxn of *M. hominis* should have Cp values between cycles 18 and 37.

8.3. Interpretation of Data

Sample 640 <i>M. hominis</i>	Sample 660 Crtl. Reaction	Channel 640 PositiveControl	Negative Control (NTC)	Result (warnings)
no amplification	detectable	amplification	negative	Negative (not detectable)
Cp < 39⁺	not relevant	amplification	negative	Positive for <i>M. hominis</i>
no amplification	not detectable	amplification	not relevant	PCR failure, repeat experiment
not relevant	not relevant	no amplification	not relevant	PCR failure, repeat experiment
not relevant	not relevant	not relevant	positive	Contamination, repeat experiment

Table 5. Typical analysis results (LightCycler® 480 II Instrument, Roche Master: Fast Start)

⁺ The cut off value is a recommendation only - this value has to be defined by the user. Compare with the lot specific values for 10 copies as reported in the Certificate of Analysis (CoA). 1-2 cycles later than observed for 10 copies corresponds to ~5 copies.

8.4. Sample Data – Typical Results




Fig.2. LightCycler® 480 II sample data for the *M. hominis* detection system.

Upper panels: Left panel Filter Combination 498-640 quantification mode (Second Derivative Maximum) with amplification curves for *M. hominis*. Right panel Filter Combination 498-640 melting analysis for *M. hominis*.

Central panels: Left panel Filter Combination 498-660 quantification mode (Second Derivative Maximum) for the control reaction. Right panel Filter Combination 498-660 melting analysis for the control reaction (not relevant for detection).

Lower panels: Left panel Filter 498-640 quantification mode (Second Derivative Maximum) for *M. hominis* and *M. genitalium*. Right panel melting analysis - identification of the species if combined with LightMix® Kit 40-0169 *M. genitalium*.

Note: Fluorescence levels may vary also depending on instrument settings. The characteristics of the curve (shape and trend of fluorescence levels) should be similar to the curve in the manual. The fluorescent curves over cycles (quantification mode) must be smooth and not zigzag shaped. The Cp values will vary from instrument to instrument by up to 2 cycles.

9. References

¹ Development of real-time PCR for detection of *Mycoplasma hominis*. Baczynska et al., 2004

10. Material Safety Data

According to OSHA 29CFR1910.1200, Commonwealth of Australia [NOHSC:1005, 1008 (1999)] and the European Union Directives 67/548/EC and 1999/45/EC any products which not contain more than 1% of a component classified as hazardous or classified as carcinogenic do not require a Material Safety Data Sheet (MSDS).

Product is not hazardous, not toxic, not IATA-restricted. Product is not from human, animal or plant origin. Product contains synthetic oligonucleotide primers and probes.

11. Version History

Events requiring changes in procedures red, mod. sequences blue

V120214	Release version
V121005	Data for the use in combination with LightMix Kit 40-0169 <i>M. genitalium</i>
V130308	Update section 6.3 pipetting volumes. cobas z 480 Instruments included. Use in combination with LightMix Kit 40-0153-32 <i>Ureaplasma urealyticum</i> . Chapters 7.3 and 7.4 as well as 8.3 and 8.4 are exchanged.
V130813	Editorial changes
V140630	Combination with kit 40-0153 removed from instructions. Alternative kit 40-0460 including <i>U. urealyticum</i> added to the manual. Primer for <i>M.hom</i> adapted to detect previously missed variants. Some isolates of <i>M.hominis</i> may exhibit a Tm of about 56°C
V150404	Change Internal Control (IC) to Extraction Control (EC)
V150505	Change to universal ¹²⁵ ECT target containing Lambda and PhHV DNA

Roche SAP order n° 06295070001

Notice to Purchaser

LightCycler[®] hybridization probes and kits produced under license agreements with Roche. The purchase price of this product includes a limited, nontransferable license under US patent claims and corresponding patent claims outside the United States, licensed from BioFire Defense, to use only this amount of the product for HybProbe assays and related processes described in said patents solely for the research and development activities of the purchaser. Diagnostic uses require a separate license from Roche.

These reagents were developed and manufactured by TIB MOLBIOL GmbH, Berlin, Germany.

