

LightMix[®] Kit *human Adenovirus* (hAdV)

Cat.-No. 40-0303-32

2015: 32rxn/vial and IC changed to Extraction Control

Kit with reagents for the detection of *human Adenovirus* DNA using the Roche Diagnostics LightCycler[®] 1.x / 2.0 / 480 II Instruments or cobas z 480 Analyzer.

Lyophilized mix of primers and probes for a total of 96 reactions with a final volume of 20 µl each.
Shipping at ambient temperature. Store protected from light at 4°C-25°C - do not freeze !

Instructions for use with the LightCycler[®] 1.x / 2.0 Instruments see pages 4-5

Instructions for use with the LightCycler[®] 480 II Instrument / cobas z 480 Analyzer see pages 6-7

1. Introduction

Adenoviruses have a linear double stranded ~ 40 kb DNA genome, containing 30-40 genes. The virus affects membranes of the respiratory tract, eyes, intestines, and urinary tract. Symptoms of respiratory illness range from the common cold to pneumonia and bronchitis. Adenoviruses account for about 10% of acute respiratory infections in children and are a frequent cause of diarrhea. The disease is usually mild and requires at most a symptomatic treatment. Immune-compromised patients are vulnerable to severe complications of Adenovirus infections. Transmission is airborne or fecal route.

The 51 serotypes of human adenovirus (HAdV) are divided into six groups named from A to F. Respiratory diseases are mainly due to virus from groups B and C, conjunctivitis to groups C and D, while the causative virus for gastroenteritis is probably limited to group F viruses (types 40 and 41 only).

| Group | Type |
|-------|--|
| A | 12, 18, 31 |
| B | 3, 7, 11, 14, 21, 34, 35 |
| C | 1, 2, 5, 6 |
| D | 8-10, 13, 15, 17, 19, 20, 22-30, 32, 33, 36-39 |
| E | 4 |
| F | 40, 41 |

Routine diagnosis is accomplished by antigen testing, serology, or PCR. Clinical samples can be swabs, feces, rectal swabs and tissue as well as blood, serum, plasma, or liquor.

2. Description

This kit provides a fast and accurate system to detect *Adenovirus* (AdV) genomic DNA from a nucleic acid extract; the kit includes a spiked extraction control working also as amplification control reaction. The kit has been verified to work with AdV types 1-51 and type 53. For a determination of the viral load in liquid samples please see section 10. Conversion Factor.

A 129 bp long fragment of the Hexon gene is amplified with several specific primers and detected with probes labeled with LightCycler[®] Red 640 (channel 640). The presence of the PCR product can be verified by running a melting curve, but without the capability of serotyping. If typing is required, the PCR product is ready to be analyzed using a Chipron LCD chip for serotype group determination.

The control reaction is based on a 278 bp long PCR fragment generated from Lambda DNA, detected with LightCycler[®] Red 690 labeled hybridization probes (channel 705). The control reaction is supplied in a separate vial and may be optionally run in a separate reaction instead of duplex PCR in order to achieve a higher sensitivity, in particular as the AdV PCR is based on a complex mixture of primers.

The former internal control (IC) has been changed to a spiked extraction control (sEC) to monitor successful extraction and demonstrate the ability to run an amplification reaction (no PCR inhibition). We recommend to use the 'Extraction Control' procedure; for compatibility with the former procedure the use as IC is also described. The extraction control target ⁿECT contains Lambda and PhHV DNA and may be used also for other LightMix kits, without adding the ECT provided with the other kit(s).

The use of a color compensation file generated with the ColorCompensation kit 40-0318 is a prerequisite to run is a prerequisite to detect the control reaction.

The supplied standard row (AdV group B) allows to determine the linear range and the detection limit of the device.

The kit must be used with 'LightCycler[®] FastStart DNA Master HybProbe' or "LightCycler[®] FastStart^{PLUS} DNA Master Hybridization Probes" (capillary and plate based LightCycler[®] Instruments).

3. Set Contents

- 3 Vials with **blue** cap containing premixed lyophilized primers and probes for 32 PCR rxns *HAdV*
- 3 Vials with **white** cap containing premixed primers and probes for each 32 control reactions

- 1 Standard row with 6 lyophilized plasmid standards of HAdV from 10^1 to 10^6 target equivalents / rxn
- 1 Sealing foil for the standard row

- 1 Vial with **white** cap with the universal Extraction Control Target (ⁿ**ECT**): 4.8×10^6 copies (total)
- 1 Vial with **black** cap containing the stabilizer for the No Template Control (NTC)

- 1 Certificate of Analysis (CoA) with lot-specific data

4. Additional Reagents and Items Required

| | |
|--|--|
| ColorCompensation HybProbe order n° 40-0318-00 | Roche Diagnostics Cat.-No. 05 997 704 001 |
| LightCycler® FastStart DNA Master HybProbe or LightCycler® FastStart DNA Master ^{PLUS} HybProbe | Cat.-No. 03 003 248 001 Cat.-No. 03 515 575 001 |
| High Pure PCR Template Preparation Kit | Cat.-No. 11 796 828 001 |
| High Pure Viral Nucleic Acid Kit | Cat.-No. 11 858 874 001 |
| LightCycler® Capillaries (20 µl) (LightCycler® 1.x / 2.0 Instruments) | Cat.-No. 04 929 292 001 |
| LightCycler® 480 Multiwell Plate 384, white (LightCycler® 480 Instrument) or LightCycler® 480 Multiwell Plate 96, white (LightCycler® 480 Instrument) | Cat.-No. 04 729 749 001 Cat.-No. 04 729 692 001 |

For use in LightCycler® 1.x Instruments with software version 3.5.3 read channel F2 instead of 640 and F3 instead of channel 705. We recommend to upgrade to software version 4.10 or higher.

5. Product Characteristics

PCR results (activation, 50 cycles and melting curve) are obtained within 50 minutes with the capillary based LightCycler® 1.x / 2.0 Instruments and within 80 minutes with '480' plate based Instruments.

Sensitivity

Reagents detect 10 copies/reaction positive control target DNA using FastStart DNA Master HybProbe.

Measuring range

The linear measuring range of the assay is 10^2 to 10^6 copies of target genomic DNA.

Storage and Stability

- Lyophilized reagents are stable for at least 6 months after shipment when stored protected from light at 4°C-25°C. See expiry date on the product label.
- **Do not freeze** lyophilized reagents.
- Dissolved positive controls must be stored at -20°C. Avoid multiple freeze-thaw cycles.
- Dissolved reagents are stable for at least 10 days when stored protected from light and refrigerated (4°C).
- Dissolved reagents can be long-term stored frozen at -20°C. Avoid multiple thaw-freeze cycles.

6. Experimental Protocol

Start programming before preparing the solutions. See the Instrument operator's manual for details.

6.1. Preparation of Parameter-Specific Reagents (PSR):

One reagent vial with a **blue** cap contains primers and probes to run 32 reactions of *Adenovirus*. One reagent vial with a **white** cap contains primers and probes to run 32 Control Reactions.

Check for the colored pellet, then **add 66 µl** PCR-grade water, mix (vortex) and spin down.

► Use **2 µl** reagent for a 20 µl PCR reaction.

| This solution is stable at least ten days when stored refrigerated at 4°C. Avoid prolonged exposure to light.

6.2. Preparation of Extraction Control Target (ECT):

Add **1,200 µl** PCR-grade water to the vial with **white** cap containing the Extraction Control Target

6.3. Sample Preparation:

Before extraction add **10 µl** of ECT (vial **white** cap) to the sample. **Skip if IC procedure is used.** Perform nucleic acid preparation as described in the protocol of the extraction kit used (see 4.1).

6.4. Preparation of No Template Extraction (NTC):

Add **900 µl** PCR-grade water to the **NTC** (vial **black** cap) and add **100 µl** of ECT (vial **white** cap).

► Use **5 µl** No Template Control DNA for a 20 µl PCR reaction.

6.5. Preparation of the Standard Row

The target DNA is provided in 6 different quantities to yield from 10^1 to 10^6 target molecules in 5 µl once resolved. Start with the lowest concentrated standard (first tube from the extended lip). Use the pipette tip to punch a hole through the sealing foil. Add **40 µl** of **NTC** (vial **black** cap) to each well. **Use water if the IC procedure is selected.** Mix the target DNA by pipetting the solution up and down 10 times.



► Use **5 µl** standard for a 20 µl PCR reaction.

| This standard solution is not long-term stable. Use only fresh prepared solutions only. After adding the target DNA to the reaction mix, use the provided sealing foil to close the vials in order to avoid contaminations.

6.6. Preparation of the Reaction Mix

In a cooled reaction tube, prepare the reaction mix by multiplying each volume for a single reaction by the number of reactions to be cycled plus one additional reaction.

| For use with the Roche FastStart ^{PLUS} Master | | For use with the Roche FastStart Master | | |
|---|--------------------|---|--------------------|--------------------|
| Single reaction IC | Single reaction EC | Component | Single reaction EC | Single reaction IC |
| 6.5 µl | 7.0 µl | water, PCR-grade (colorless cap, provided with the Roche Master kit) | 7.4 µl | 6.9 µl |
| -- | -- | Mg ²⁺ solution 25 mM (blue cap, provided with the Roche FastStart kit) | 1.6 µl | 1.6 µl |
| 2.0 µl | 2.0 µl | PSR mix (parameter specific reagents, see 6.1.) | 2.0 µl | 2.0 µl |
| 2.0 µl | 2.0 µl | Primers and probe mix for the IC/EC | 2.0 µl | 2.0 µl |
| 0.5 µl | -- | ECT (white cap, DNA control target, see 6.2.) | -- | 0.5 µl |
| 4.0 µl | 4.0 µl | Roche Master (red cap, for preparation see Roche manual) | 2.0 µl | 2.0 µl |

15.0 µl

Volume of reaction mix

15.0 µl

Table 1

Mix gently, spin down and **transfer 15 µl** of the reaction mix to a capillary or well.

Add 5 µl of sample or standard to each capillary or well for a final reaction volume of 20 µl. Close the capillaries / attach a foil to the multiwell plate and seal, and spin down.

Start run.

7. LightCycler® 1.x / 2.0 Instruments

7.1. Programming

The protocol consists of four program steps

- 1: Denaturation: sample denaturation and enzyme activation
- 2: Cycling: PCR-amplification of the target DNA
- 3: Melting: melting curve analysis for identification of the PCR product derived from the target DNA
- 4: Cooling: cooling the instrument

| Program Step: | Denaturation | Cycling | | | Melting | | | Cooling |
|---------------------|--------------|---------------------|----------|----------|---------------------|----------|----------|----------|
| Parameter | | | | | | | | |
| Analysis Mode | None | Quantification mode | | | Melting Curves mode | | | None |
| Cycles | 1 | 50 | | | 1 | | | 1 |
| Target [°C] | 95 | 95 | 62 | 72 | 95 | 40 | 85 | 40 |
| Hold [hh:mm:ss] | 00:10:00 | 00:00:05 | 00:00:05 | 00:00:15 | 00:00:20 | 00:00:20 | 00:00:00 | 00:00:30 |
| Ramp Rate [°C/s] | 20 | 20 | 20 | 20 | 20 | 20 | 0.2 | 20 |
| Sec Target [°C] | - | - | 55 | - | - | - | - | - |
| Step Size [°C] | - | - | 0.5 | - | - | - | - | - |
| Step Delay (Cycles) | - | - | 1 | - | - | - | - | - |
| Acquisition Mode | None | None | Single | None | None | None | Cont | None |

Table 2

7.2. Data Analysis

Switch the color compensation mode on. If this mode is not enabled run the color compensation program. Follow the instructions in the manual Color Compensation Kit HybProbe 530/640/690.

Perform data analysis, as described in the LightCycler® Instrument operator's manual.

We recommend using the Second Derivative Maximum method (Automated (F'' max)). The cycle number of the Crossing Point (Cp) of each sample is calculated automatically. The Fit Points method is more-error prone due to the user's influence.

View *Adenovirus* data in channel 640, Quantification mode (LightCycler® 2.0 Instrument). The negative control (NTC) must show no signal.

For the identification of the PCR product view *Adenovirus* data in channel 640 Melting Curves mode.

For the Control Reaction view channel 705 data. The negative control and the low-concentrated *Adenovirus* DNA samples (10 to 1,000 copies) should show an amplification curve for the EC/IC with a Cp at approximately cycle 25-31.

The provided standard row of cloned DNA with concentrations in the range from 10⁶ copies/rxn to 10¹ copies/rxn of *Adenovirus* should have Cp values between cycles 18 and 36.

7.3. Interpretation of Data

| Sample 640 <i>Adenovirus</i> | Sample 705 Ctrl Reaction | Channel 640 PositiveControl | Negative Control (NTC) | Result (warnings) |
|---------------------------------|-----------------------------|--------------------------------|---------------------------|--|
| no amplification | detectable | amplification | negative | Negative (not detectable) |
| Cp < 37⁺ | not relevant | amplification | negative | Positive for <i>Adenovirus</i> |
| no amplification | not detectable | amplification | not relevant | PCR failure , repeat experiment |
| not relevant | not relevant | not relevant | positive | Contamination , repeat experiment |

Table 3. Typical analysis results (LightCycler® 2.0 Instrument, Roche Master: FastStart)

⁺ The cut off value is a recommendation only - this value has to be defined by the user. Compare with the lot specific values for 10 copies as reported in the Certificate of Analysis (CoA). 1-2 cycles later than observed for 10 copies corresponds to ~5 copies

7.4. Sample Data – Typical Results




Fig.1. LightCycler® 2.0 sample data for the Adenovirus detection system.

Upper panels: Left panel channel 640 quantification mode (Second Derivative Maximum) with amplification curve for AdV.

Middle panels: Left panel channel 640 quantification mode (Second Derivative Maximum) for the samples.

Right panel channel 640 melting peaks for Adenovirus groups A to F.

Lower panels: Left panel channel 705 quantification mode (Second Derivative Maximum) with amplification curves for the control reaction. 705 melting analysis for the control reaction (not relevant for detection).

* **Note:** Fluorescence levels depend on instrument settings and may vary. The characteristics of the curve (shape and trend of fluorescence levels) should be similar to the curve in the manual. The fluorescent curves over cycles (quantification mode) must be smooth and not zigzag shaped. The Cp values will vary from instrument to instrument by up to 2 cycles, while the distance between two dilution steps should be relative constant (delta Cp). The Cp values described in this manual (chart text) have been obtained using the supplied standard row; the entire set of curves can be horizontal shifted.

8. LightCycler® 480 II Instrument / cobas z 480 Analyzer

8.1. Programming

The protocol consists of four program steps

- 1: Denaturation: sample denaturation and enzyme activation
- 2: Cycling: PCR-amplification of the target DNA
- 3: Melting: melting curve analysis for identification of the PCR product derived from the target DNA
- 4: Cooling: cooling the instrument

Detection Format:

LightCycler® 480 II Instrument: 465-510, 498-640, 498-660

cobas z 480 Analyzer: 465-510, 498-645, 498-700

| Program Step: | Denaturation | Cycling | | | Melting | | | Cooling |
|-----------------------|--------------|---------------------|----------|----------|---------------------|----------|------------|----------|
| Parameter | | | | | | | | |
| Analysis Mode | None | Quantification mode | | | Melting Curves mode | | | None |
| Cycles | 1 | 50 | | | 1 | | | 1 |
| Target [°C] | 95 | 95 | 62 | 72 | 95 | 40 | 85 | 40 |
| Hold [hh:mm:ss] | 00:10:00 | 00:00:05 | 00:00:05 | 00:00:15 | 00:00:30 | 00:02:00 | 00:00:00 | 00:00:30 |
| Ramp Rate [°C/s] 96 | 4.4 | 4.4 | 2.2 | 4.4 | 4.4 | 1.5 | - | 1.5 |
| Ramp Rate [°C/s] 384 | 4.6 | 4.6 | 2.4 | 4.6 | 4.6 | 2.0 | - | 2.0 |
| Sec Target [°C] | - | - | 55 | - | - | - | - | - |
| Step Size [°C] | - | - | 0.5 | - | - | - | - | - |
| Step Delay (Cycles) | - | - | 1 | - | - | - | - | - |
| Acquisition Mode | None | None | Single | None | None | None | Continuous | None |
| Acquisitions [per °C] | - | - | - | - | - | - | 1 | - |

Table 4

8.2. Data Analysis

Note: cobas z 480 Analyzer signal levels are about 50% compared to LightCycler® 480 II results.

Switch the color compensation mode on. If this mode is not enabled run the color compensation program. Follow the instructions in the manual 'LightMix® Kit – Color Compensation HybProbe.

Perform data analysis, as described in the LightCycler® Instrument operator's manual.

We recommend using the Second Derivative Maximum method (Automated (F" max)). The cycle number of the Crossing Point (Cp) of each sample is calculated automatically.

View *Adenovirus* data with Filter Combination 498-640, Quantification mode. The negative control (NTC) must show no signal. For the identification of the PCR product view *Adenovirus* data with Filter Combination 498-640 Melting Curves mode.

If the control reaction is used, view data with Filter Combination 498-660 Quantification mode. The negative control and the low-concentrated *Adenovirus* DNA samples (10 to 1,000 copies) should show an amplification curve for the EC/IC with a Cp at approximately cycle 25-31.

The provided standard row of cloned and purified DNA with concentrations in the range from 10⁶ copies/rxn to 10¹ copies/rxn of *Adenovirus* should have Cp values between cycles 18 and 36.

8.3. Interpretation of Data

| Sample 640 <i>Adenovirus</i> | Sample 660 Ctrl Reaction | Channel 640 PositiveControl | Negative Control (NTC) | Result (warnings) |
|---------------------------------|-----------------------------|--------------------------------|---------------------------|----------------------------------|
| no amplification | detectable | amplification | negative | Negative (not detectable) |
| Cp < 36 ⁺ | not relevant | amplification | negative | Positive for <i>Adenovirus</i> |
| no amplification | not detectable | amplification | not relevant | PCR failure, repeat experiment |
| not relevant | not relevant | not relevant | positive | Contamination, repeat experiment |

Table 5. Typical analysis results (LightCycler® 480 Instrument, Roche Master: FastStart)

⁺ The cut off value is a recommendation only - this value has to be defined by the user. Compare with the lot specific values for 10 copies as reported in the Certificate of Analysis (CoA). 1-2 cycles later than observed for 10 copies corresponds to ~5 copies

8.3. Sample Data – Typical Results




Fig.2. LightCycler® 480 II sample data for the Adenovirus detection system.

Upper panels: Left panel Filter Combination 498-640 quantification mode (Second Derivative Maximum) with amplification curves for Adenovirus (Fast Start). Right panel Filter Combination 498-640 amplification curves for Adenovirus (Fast Start^{Plus}).

Middle panels: Left panel Filter Combination 498-640 quantification mode (Second Derivative Maximum) with amplification curves for Adenovirus (Group A-F). Right panel Filter Combination 498-640 melting peaks for Adenovirus (Group A-F).

Lower panels: Left panel Filter Combination 498-660 amplification curves (Second Derivative Maximum) for the control reaction. Right panel Filter Comb. 498-660 melting analysis and melting) for the control reaction (not relevant for detection).

* **Note:** Fluorescence levels may vary. The characteristics of the curve (shape and trend) should be similar to the curve shown. The fluorescent curves over cycles (quantification mode) must be smooth and not zigzag shaped. The Cp values may vary between instruments by up to 2 cycles, while the delta Cp between dilution steps should be constant. The Cp values described in this manual (chart text) have been obtained using the supplied standard row.

9. Conversion Factor

The amount of virus per sample (Viral Load) is usually reported in copies / ml while PCR reports copies per reaction. The conversion factor between both depends on sample and extraction volumes. Extraction starts usually from less than one milliliter and PCR does not use the total extracted volume.

The viral load (VL) can be calculated using the following general formula:

$$\text{VL [copies/ml]} = \text{MV (Measured Value)} \times \text{EVF} \times \text{SVF}$$

Example: Extracting 200 µl clinical sample results in an extraction volume factor (EVF) of **5**. Using 5 µl sample from a total elution volume of 100 µl results in a sample volume factor (SVF) of **20**, resulting in a conversion factor of **100**:

$$\text{Viral Load [copies/ml]} = \text{Measured Value [copy number per reaction]} \times 100$$

10. References (Use of this kit)

Adenovirus Disease after Kidney Transplantation: Course of Infection and Outcome in Relation to Blood Viral Load and Immune Recovery. Watcharananana et al. (2011)

Swabbing for resp. viral infections in older patients: a comparison of rayon and nylon flocked swabs. Hernes et al. (2011)

Alt Solunum Yolu Enfeksiyonu Olan Çocuklarda Adenovirusların Araştırılması. Saglic et al. Mikrobiyol Bul (2013)

11. Material Safety Data

According to OSHA 29CFR1910.1200, Australia [NOHSC:1005, 1008 (1999)] and the European Union Directives 67/548/EC and 1999/45/EC any products which not contain more than 1% of a component classified as hazardous or classified as carcinogenic do not require a MSDS.

Product is not hazardous, not toxic, not IATA-restricted. Product is not from human, animal or plant origin. Product contains synthetic oligonucleotide primers and probes.

12. Version History Red notes mark changes in procedures, blue modification of sequences

| | |
|---------|---|
| V111208 | Correction of Cp values for the IC in LC 2.0 |
| V130813 | z 480 included, Conversion Factor, manual condensed to 8 pages |
| V150505 | Kit changed from 6 vials 16 rxns to 3 vials 32 rxns Internal control changed to extraction control (EC) Kit equipped with the universal ⁿECT extraction control target |
| V160330 | Cp values control reaction corrected, Storage conditions, Section 9 wording |

Roche SAP order n° 05552451001

Notice to Purchaser

LightCycler[®] hybridization probes and kits produced under license agreements with Roche. The purchase price of this product includes a limited, nontransferable license under US patent claims and corresponding patent claims outside the United States, licensed from BioFire Defense, to use only this amount of the product for HybProbe assays and related processes described in said patents solely for the research and development activities of the purchaser. Diagnostic uses require a separate license from Roche.

These reagents were developed and manufactured by TIB MOLBIOL GmbH, Berlin, Germany.