

LightMix[®] Kit *Mycoplasma pneumoniae* EC

Cat.-No. 40-0221-32

Internal Control (IC) changed to spiked Extraction Control (sEC)

Kit with reagents for detection of *Mycoplasma pneumoniae* genomic DNA using Roche Diagnostics LightCycler[®] 1.x / 2.0 / 480 II Instruments or cobas z 480 Analyzer.

Lyophilized mix of primers and probes for a total of 96 reactions with a final volume of 20 µl each.

Store protected from light at room temperature (18-25°C), do NOT freeze!

Instructions for use with the LightCycler[®] 1.x / 2.0 Instruments see pages 4-5

Instructions for use with the LightCycler[®] 480 II / cobas z 480 Instruments see pages 6-7

1. Introduction

Mycoplasma pneumoniae, like the other members of the genus, lacks a bacterial cell wall and has one of the smallest bacterial genomes known. For several essential biomolecules such as fatty acids, amino acids, precursors for nucleic acids and cholesterol *M. pneumoniae* depend on the host.

M. pneumoniae can cause tracheobronchitis and primary atypical pneumonia with clinical symptoms such as nonproductive cough sore throat, low-grade fever, and middle-ear involvement. Together with the bacteria *Chlamydia pneumoniae* and *Legionella spp./pneumophila* it belongs to the common agents causing 'Atypical Pneumonia'.

Due to the lack of a cell wall all B-lactam antibiotics (e.g. penicillin, cycloserine) are not effective. Other antibiotic therapies which are usually ineffective on bacteria such as the use of polyenes acting on the cholesterol of the cell membrane allow for an avenue to fight *M. pneumoniae*. However, the most efficient way for a host to fight off a *Mycoplasma* infection is through its immune response.

Common targets for PCR based detection of *M. pneumoniae* are the 16S RNA gene, the ATPase operon orf521, the repetitive element repMp1, or the cytoadhesin protein gene P1^{1,2}.

¹ Evaluation of Chlamydia pneumoniae and Mycoplasma pneumoniae as Etiologic Agents of Persistent Cough in Adolescents and Adults. Wadowsky et al., JCM (2002)

² Development and evaluation of a loop-mediated isothermal amplification assay for rapid detection of Mycoplasma pneumoniae. Saito et al., JMM (2005)

2. Description

This kit provides a fast and accurate system to detect and identify *Mycoplasma pneumoniae* genomic DNA in a nucleic acid extract. A 114 bp fragment of the cytoadhesin protein P1 gene is amplified with specific primers and detected with probes labeled with LightCycler[®] Red 640 (detected in channel 640).

The control reaction generates an additional product of 247 bp, using hybridization probes labeled with LightCycler[®] Red 690 (recorded in channel 705). This second PCR has no visible impact on the Mycoplasma-specific reaction and will even fail in the presence of higher amounts of target (1,000 copies and more) while displaying an amplification signal in negative and low-concentrated samples.

The former internal control (IC) has been changed to a spiked extraction control (sEC) to monitor a successful extraction and demonstrate the ability to run a PCR reaction (absence of inhibition).

We recommend to use the 'Extraction Control' procedure; for compatibility with the former procedure the use as IC is also described. Target and control primer/probe sequences remained unchanged.

Note: With a MagNA Pure Compact extraction the recovery rate for the EC target can be very low or the DNA target even gets lost.

The use of a color compensation file generated with the LightMix[®] Color Compensation HybProbe 530/640/690 kit (order no. 40-0318-00) is a prerequisite to detect the control reaction.

The supplied standard row DNA allows to determine the linear range of the reaction and to estimate the quantity of the target sequence in unknown samples.

The kit is tested with 'LightCycler[®] FastStart DNA Master HybProbe' only.

3. Set Contents

- 3 Vials with **green** cap containing premixed primers and probes for 32 PCR reactions
- 3 Vials with **white** cap containing premixed primers and probes for 32 control reactions
- 1 Vial with white cap containing Extraction Control Target (ECT): 4.8×10^6 copies (total)
- 1 Vial with black cap containing the stabilizer for preparation of the 'No Template Control' (NTC)
- 1 Standard row with 6 lyophilized standards *M. pneumoniae* $10^1 - 10^6$ target equivalents / rxn
- 1 Sealing foil for the standard row
- 1 Certificate of Analysis (CoA) with lot-specific data

4. Additional Reagents and Items Required

Color Compensation HybProbe order n° 40-0318-00	Roche Diagnostics Cat.-No. 05 997 704 001
LightCycler® FastStart DNA Master HybProbe	Cat.-No. 03 003 248 001
LightCycler® Capillaries (20 µl) (LightCycler® 1.x / 2.0 Instruments)	Cat.-No. 04 929 292 001
LightCycler® 480 Multiwell Plate 384, white (LightCycler® 480 Instrument) or LightCycler® 480 Multiwell Plate 96, white (LightCycler® 480 Instrument)	Cat.-No. 04 729 749 001 Cat.-No. 04 729 692 001

For use in LightCycler® 1.x Instruments with software version 3.5.3 read channel F2 instead of channel 640, channel F3 instead of channel 705 and channel F1 instead of channel 530 for detection. We recommend upgrading LightCycler® 1.x Instruments to software version 4.1.

4.1. Optional Additional Reagents

High Pure PCR Template Preparation Kit	Cat.-No. 11 796 828 001
--	-------------------------

5. Product Characteristics

PCR results are obtained within 50 minutes (50 cycles and melting curve) with the LightCycler® 1.x / 2.0 Instruments and within 80 minutes (50 cycles and melting curve) with plate based Instruments.

Sensitivity

These reagents detect 10 copies of *M. pneumoniae* DNA per reaction, using Roche 'LightCycler® FastStart DNA Master HybProbe' (in an exemplary system, using cloned targets as reference).

Measuring range

The linear measuring range of the assay is 10^2 to 10^6 copies / reaction of *M. pneumoniae* DNA using the Roche 'LightCycler® FastStart DNA Master HybProbe'.

Storage and Stability

- Lyophilized reagents are stable for at least 6 months after shipment when stored protected from light at room temperature (18-25°C). See expiry date on the product label.
- **Do not freeze** lyophilized reagents.
- Dissolved reagents are stable for at least 10 days when stored protected from light and refrigerated (4°C).

6. Experimental Protocol

Start programming before preparing the solutions. See the instrument operator's manual for details.

6.1. Preparation of Parameter-Specific Reagents (PSR):

One reagent vial with a **green** cap contains primers and probes to run 32 reactions *M. pneumoniae*. One reagent vial with a **white** cap contains primers and probes to run 32 reactions control reaction.

Check for the colored pellet, then **add 66 µl** PCR-grade water, mix (vortex) and spin down.

► **Use 2 µl** reagent for a 20 µl PCR reaction.

| This solution is stable at least ten days when stored refrigerated at 4°C. Avoid prolonged exposure to light.

6.2. Preparation of Extraction Control Target (ECT):

Add **1,200 µl** PCR-grade water to the vial with **white** cap containing the Extraction Control Target

6.3. Sample Preparation:

Before extraction add **10 µl** of ECT (vial **white** cap) to the sample. **Skip if IC procedure is used.** Perform nucleic acid preparation as described in the protocol of the extraction kit used (see 4.1).


6.4. Preparation of No Template Extraction (NTC):

Add **900 µl** PCR-grade water to the **NTC** (vial **black** cap) and add **100 µl** of ECT (vial **white** cap).

► **Use 5 µl** No Template Control DNA for a 20 µl PCR reaction.

6.5. Preparation of the Standard Row:

The target DNA is provided in 6 different quantities to yield from 10^1 to 10^6 target molecules in 5 µl once resolved. Start with the lowest concentrated standard (first tube from the extended lip). Use the pipette tip to punch a hole through the sealing foil. Add **40 µl** of **NTC** (vial **black** cap) to each well. **Use water if the IC procedure is selected.** Mix the target DNA by pipetting the solution up and down 10 times.



► **Use 5 µl** standard DNA for a 20 µl PCR reaction.

| This standard solution is not long-term stable and will lose sensitivity under prolonged storage. After adding the target DNA to reaction mix, use the provided sealing foil to close the vials in order to avoid contaminations.

6.6. Preparation of the Reaction Mix

Include Positive Controls and at least one 'No Template Control' (NTC). In a cooled tube, prepare the reaction mix by multiplying the single reaction volumes by the number of reactions plus one reserve :

sEC Procedure	For use with the Roche FastStart Master	IC Procedure
Single reaction	Component	Single reaction
7.4 µl	Water, PCR-grade (colorless cap, provided with the Roche Master kit)	6.9 µl
1.6 µl	Mg ²⁺ solution 25 mM (blue cap, provided with the Roche FastStart kit)	1.6 µl
2.0 µl	PSR mix (parameter specific reagents, primers and probes, see 6.1.)	2.0 µl
2.0 µl	Control Reaction [see 6.1]	2.0 µl
---- µl	ECT Control Target (vials white cap)	0.5 µl
2.0 µl	Roche Master (red cap, for preparation see Roche manual)	2.0 µl
15.0 µl	Volume of reaction mix	15.0 µl

Table 1

To run the assay without the control reaction substitute ECT with 0.5 µl PCR-grade water.

Mix gently, spin down and **transfer 15 µl** of the reaction mix to a capillary or well.

Add 5 µl of sample or standard to each capillary or well for a final reaction volume of 20 µl. Close the capillaries / attach a foil to the multiwell plate and seal, and spin down.

Start run.

7. LightCycler® 1.x / 2.0 Instruments

7.1. Programming

The protocol consists of four program steps

- 1: Denaturation: sample denaturation and enzyme activation
- 2: Cycling: PCR-amplification of the target DNA
- 3: Melting: melting curve analysis for identification of the PCR product derived from the target DNA
- 4: Cooling: cooling the instrument

Program Step:	Denaturation	Cycling			Melting			Cooling
Parameter								
Analysis Mode	None	Quantification mode			Melting Curves mode			None
Cycles	1	50			1			1
Target [°C]	95	95	62	72	95	40	85	40
Hold [hh:mm:ss]	00:10:00	00:00:05	00:00:05	00:00:15	00:00:20	00:00:20	00:00:00	00:00:30
Ramp Rate [°C/s]	20	20	20	20	20	20	0.2	20
Sec Target [°C]	-	-	55	-	-	-	-	-
Step Size [°C]	-	-	0.5	-	-	-	-	-
Step Delay (Cycles)	-	-	1	-	-	-	-	-
Acquisition Mode	None	None	Single	None	None	None	Cont	None

Table 2

(Melting not relevant for detection)

7.2. Data Analysis

Switch the color compensation mode on. If this mode is not enabled run the color compensation program. Follow the instructions in the manual Color Compensation Kit HybProbe 530/640/690.

Perform data analysis, as described in the Instrument operator's manual.

We recommend using the Second Derivative Maximum method (Automated (F'' max)). The cycle number of the Crossing Point (Cp) of each sample is calculated automatically. The Fit Points method is more-error prone due to the user's influence.

View *M. pneumoniae* data in channel 640, Quantification mode. The negative control (NTC) must show no signal. For the identification of the PCR product view *M. pneumoniae* data in channel 640 Melting Curves mode.

For the Control Reaction view channel 705 data. The negative control and the low-concentrated DNA samples (10 to 1,000 copies) will show an amplification curve with a Cp at approximately cycle 25-28.

The provided standard row of cloned *M. pneumoniae* target DNA with concentrations from 10⁶ to 10 copies/reaction target should yield Cp values between cycles 16 and 34 (see figure 1).

For use in LightCycler® 1.x / 2.0 Instruments select channel F2 instead of channel 640, channel F3 instead of channel 705 for detection.

7.3. Sample Data – Typical results




Fig.1. LightCycler® 2.0 sample data for the *Mycoplasma pneumoniae* detection system.

Upper panels: Left panel channel 640 quantification mode (Second Derivative Maximum) with amplification curves for *Mycoplasma pneumoniae*. Right panel channel 640 melting analysis for *Mycoplasma pneumoniae* (not relevant for detection).
Lower panels: Left: channel 705 quantification mode (Sec. Der. Maximum). Right: channel 705 melting peaks control reaction.

* **Note:** Fluorescence levels depend on instrument settings and may vary. The characteristics of the curve (shape and trend of fluorescence levels) should be similar to the curve in the manual. The fluorescent curves over cycles (quantification mode) must be smooth and not zigzag shaped. The Cp values will vary from instrument to instrument by up to 2 cycles, while the distance between two dilution steps should be relative constant (delta Cp). The Cp values described in this manual (chart text) have been obtained using the supplied standard row; the entire set of curves can be horizontal shifted.

7.4. Interpretation of Data

Sample 640 <i>M. pneumoniae</i>	Sample 705 Control Reaction	Channel 640 PositiveControl	Negative Control (NTC)	Result (warnings)
no amplification	detectable	amplification	negative	Negative (not detectable)
Cp < 36 ⁺	not relevant	amplification	negative	Positive for <i>M. pneumoniae</i>
no amplification	not detectable	amplification	not relevant	PCR failure, repeat experiment
not relevant	not relevant	no amplification	not relevant	PCR failure, repeat experiment
not relevant	not relevant	not relevant	positive	Contamination, repeat experiment

Table 3. Typical analysis results (LightCycler® 2.0 Instrument, Roche Master: FastStart)

⁺ The cut off value is a recommendation only - this value has to be defined by the user. Compare with the lot specific values for 10 copies as reported in the Certificate of Analysis (CoA). 1-2 cycles later than observed for 10 copies corresponds to ~5 copies.

8. LightCycler® 480 II Instrument and cobas z 480 Analyzer

8.1. Programming

The protocol consists of four program steps

- 1: Denaturation: sample denaturation and enzyme activation
- 2: Cycling: PCR-amplification of the target DNA
- 3: Melting: melting curve analysis for identification of the PCR product derived from the target DNA
- 4: Cooling: cooling the instrument

Detection Format:

LightCycler® 480 II Instrument: 465-510, 498-640, 498-660

cobas z 480 Analyzer: 465-510, 498-645, 498-700

Program Step:	Denaturation	Cycling			Melting			Cooling
Parameter								
Analysis Mode	None	Quantification mode			Melting Curves mode			None
Cycles	1	50			1			1
Target [°C]	95	95	62	72	95	40	85	40
Hold [hh:mm:ss]	00:10:00	00:00:05	00:00:05	00:00:15	00:00:30	00:02:00	00:00:00	00:00:30
Ramp Rate [°C/s] 96	4.4	4.4	2.2	4.4	4.4	1.5	-	1.5
Ramp Rate [°C/s] 384	4.6	4.6	2.4	4.6	4.6	2.0	-	2.0
Sec Target [°C]	-	-	55	-	-	-	-	-
Step Size [°C]	-	-	0.5	-	-	-	-	-
Step Delay (Cycles)	-	-	1	-	-	-	-	-
Acquisition Mode	None	None	Single	None	None	None	Continuous	None
Acquisitions [per °C]	-	-	-	-	-	-	1	-

Table 4

(Melting not relevant for detection)

8.2. Data Analysis

Note: cobas z 480 Analyzer signal levels are about 50% compared to LightCycler® 480 II results.

Switch the color compensation mode on. If this mode is not enabled run the color compensation program. Follow the instructions in the manual of the Color Compensation HybProbe kit.

Perform data analysis, as described in the Instrument operator's manual.

We recommend using the Second Derivative Maximum method (Automated (F" max). The cycle number of the Crossing Point (Cp) of each sample is calculated automatically. The Fit Points method is more-error prone due to the user's influence.

View *M. pneumoniae* data with Filter Combination 498-640 Quantification mode. The negative control (NTC) must show no signal.

For the Control Reaction view Filter Combination 498-660 data, Quantification mode. The negative control and the low-concentrated DNA samples (10 to 1,000 copies) will show an amplification curve with a Cp at approximately cycle 26-30.

The provided standard row of cloned *M. pneumoniae* target DNA with concentrations from 10⁶ to 10 copies/reaction target should yield Cp values between cycles 17 and 35 (see figure 2).

8.3. Sample Data – Typical Results




Fig.2. LightCycler® 480 II sample data for the *Mycoplasma pneumoniae* detection system.

Upper panels: Left panel Filter Combination 498-640 quantification mode (Second Derivative Maximum) with amplification curves for *M. pneumoniae*. Right panel Filter Combination 498-640 melting analysis for *M. pneumoniae* (not relevant for detection).

Lower panels: Control Reaction. Left: Filter 498-660 quantification mode (Sec. Derivative Maximum). Right : melting analysis.

* **Note:** Fluorescence levels depend on instrument settings and may vary. The characteristics of the curve (shape and trend of fluorescence levels) should be similar to the curve in the manual. The fluorescent curves over cycles (quantification mode) must be smooth and not zigzag shaped. The Cp values will vary from instrument to instrument by up to 2 cycles, while the distance between two dilution steps should be relative constant (delta Cp). The Cp values described in this manual (chart text) have been obtained using the supplied standard row; the entire set of curves can be horizontal shifted.

8.4. Interpretation of Data

Sample 640 <i>M. pneumoniae</i>	Sample 660 Control Reaction	Channel 640 Pos.Control	Negative Control (NTC)	Result (warnings)
no amplification	detectable	amplification	negative	Negative (not detectable)
Cp < 37	not relevant	amplification	negative	Positive for <i>M. pneumoniae</i>
no amplification	not detectable	amplification	not relevant	PCR failure, repeat experiment
not relevant	not relevant	no amplification	not relevant	PCR failure, repeat experiment
not relevant	not relevant	not relevant	positive	Contamination, repeat experiment

Table 5. Typical analysis results (LightCycler® 480 II Instrument, Roche Master: FastStart)

* The cut off value is a recommendation only - this value has to be defined by the user. Compare with the lot specific values for 10 copies as reported in the Certificate of Analysis (CoA). 1-2 cycles later than observed for 10 copies corresponds to ~5 copies.

10. Material Safety Data

According to OSHA 29CFR1910.1200, Commonwealth of Australia [NOHSC:1005, 1008 (1999)] and the European Union Directives 67/548/EC and 1999/45/EC any products which not contain more than 1% of a component classified as hazardous or classified as carcinogenic do not require a Material Safety Data Sheet (MSDS).

Product is not hazardous, not toxic, not IATA-restricted. Product is not from human, animal or plant origin. Product contains synthetic oligonucleotide primers and probes.

11. Version History Red notes mark events require changed procedures, blue mod. sequences

V100818	Editorial changes
V120201	Different IC
V130308	Improved reverse primer and change back to IC as in V100818 The test shows earlier Cp values and a better sensitivity. Change to 32 rxns per vial Chapter 4: Roche Color Compensation kit reference removed Chapter 8: Instrument cobas z 480 included. References and MSDS included
V130813	Editorial changes, cut-off values (recommendations)
V140909	IC changed to spiked Extraction Control
V150101	MagNA pure Compact may fail to recover the sEC extraction target.

Roche SAP order n° 05997780001



Notice to Purchaser

LightCycler® hybridization probes, Research-use and diagnostic-use kits are produced under license from Roche. The purchase price of this product includes limited, nontransferable rights under the running-royalty component to use only this amount of the product to practice the Polymerase Chain Reaction ("PCR") and related processes. No right to perform or offer commercial services of any kind using PCR, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is hereby granted by implication or estoppel.

These reagents were developed and manufactured by TIB MOLBIOL GmbH, Berlin, Germany.