

## LightMix® 480HT Scrapie Susceptibility Mutation Detection

Cat.-No. 40-0297-96

Tricolor LightMix® for the detection of ovine prion protein (PRNP) gene mutations in the codons 136, 141, 154 and 171 in a single reaction using the LightCycler® 480 Instrument.

Lyophilized mix of primers and probes for a total of 96 reactions with a final volume of 20 µl each – **store protected from light at room temperature (18-25°C), do NOT freeze!**

### 1. Set contents

- 1 Vial containing premixed and lyophilized primers and hybridization probes for 96 reactions

### 2. Experimental protocol

#### (A) Preparation of parameter-specific reagents (96 reactions):

One reagent vial contains all primers and probes to run 96 LightCycler® reactions.

Add 100 µl PCR-grade water to the reagent vial, mix the solution (vortex) and spin down.

► Use 1 µl **reagent** for a 20 µl PCR reaction.

| This solution is stable for three days or longer if stored refrigerated at 4°C. Avoid prolonged exposure to light.

#### (B) Preparation of the LightCycler® reaction mix

For use with the Roche FastStart kit		
Single reaction	Component	final
14.2 µl	water, PCR-grade (colorless cap, provided with the Roche FastStart or FastStart <sup>PLUS</sup> kit)	--
2.0 µl	FastStart mix (vial 1 (red cap), combined from vials 1a and 1b, see Roche manual)	1x
0.8 µl	Mg <sup>2+</sup> solution 25 mM (blue cap, provided with the Roche FastStart kit)	2 mM
1.0 µl	<b>reagent</b> mix (parameter specific reagents containing primers and probes, see <b>A</b> )	--
2.0 µl	sample or reference	--
20.0 µl	final volume	

Mix gently, spin down and transfer to a LightCycler® 480 Multiwell Plate.

### 3. Programming (LightCycler® 480 Instrument)

**Block Type:** 96 or 384      **Detection Format:** 450 – 500  
 483 – 533  
 450 – 610  
 450 – 640  
 483 – 670

Program:	Denaturation	Cycling			Melting			Cooling
<b>Parameter</b>		<b>Quantification</b>			<b>Melting Curves</b>			<b>None</b>
Analysis Mode	None							None
Cycles	1	45			1			1
Segment	1	1	2	3	1	2	3	1
Target [°C]	95	95	60	72	95	40	75	40
Hold [hh:mm:ss]	00:10:00	00:00:10	00:00:10	00:00:15	00:00:30	00:02:00	00:00:00	00:00:30
Ramp Rate [°C/s] <b>96</b>	4.4	4.4	2.2	4.4	4.4	1.5	-	1.5
Ramp Rate [°C/s] <b>384</b>	4.6	4.6	2.4	4.6	4.6	2.0	-	2.0
Acquisition Mode	None	None	Single	None	None	None	Continu.	None
Acquisitions [per °C]							1	


## 4. Data analysis

For data analysis use the melting curve analysis section.

Select Filter Combination 483-670 for the analysis of mutations in the codon 136 and 141 (**no color compensation necessary**).

Select Filter Combination 450-640 for the analysis of mutations in the codon 154 and Filter Combination 450-610 for the analysis of mutations in the codon 171. Switch the color compensation mode on. Channels to compensate: 450-500, 483-533, 450-610, 450-640, 483-670.

## 5. Sample data - typical results



**Fig.1. Sample data for the *ov* PRNP genotype detection system**

**Note:** Observed temperatures can differ due to buffer components (salt) coming from the extraction process and possibly due to different instrument settings. However, the temperature differences between the individual types are rather constant. Use known (pretyped) samples for calibration.

#### **Notice to Purchaser**

A license under U.S. Patents 4,683,202, 4,683,195 and 4,965,188 or their foreign counterparts, owned by Hoffmann-La Roche Inc. and F. Hoffmann-La Roche Ltd ("Roche"), has an up-front fee component and a running-royalty component. The purchase price of this product includes limited, nontransferable rights under the running-royalty component to use only this amount of the product to practice the Polymerase Chain Reaction ("PCR") and related processes described in said patents solely for the research and development activities of the purchaser when this product is used in conjunction with a thermal cycler whose use is covered by the up-front fee component. Rights to the up-front fee component must be obtained by the end user in order to have a complete license. These rights under the upfront fee component may be purchased from Perkin-Elmer or obtained by purchasing an authorized thermal cycler. No right to perform or offer commercial services of any kind using PCR, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is hereby granted by implication or estoppel. Further information on purchasing licenses to practice the PCR process for research applications may be obtained by contacting the Director of Licensing at The Perkin-Elmer Corporation, 850 Lincoln Center Drive, Foster City, California 94404 or at Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501. The purchase of this product does not convey any right for its use in clinical diagnostic applications. No rights for TaqMan technology under U.S. Patents 5,210,015 and 5,487,972 are hereby conveyed.

Reverse primer was changed December 2006 based on the report of a deviation due to a SNP at codon 176; Version 317.

These reagents were developed and manufactured by TIB MOLBIOL GmbH, Berlin, Germany.  
LightCycler<sup>®</sup> hybridization probes produced under license from Roche Diagnostics.

