



Instructions for life science research use only. Not tested for use in diagnostic procedures. For *in vitro* use only.



Instructions For Use

LightMix[®] Modular PhHV Internal Control

660

Cat.-No. 66-0625-96

Roche SAP n° 07 093 853 001

Kit with reagents for 96 PCR reactions 20 µl for detection of PhHV [lyophilized]

1. Content, Storage and Expiry

1 Vial green cap 96 reactions IC (lyophilized)

Storage at Arrival:

Store cooled or at ambient temperature
Do **not** freeze the lyophilized reagents.

- Lyophilized kits are stable for at least 6 months (4°C to 25°C in the dark). See lot-specific expiry date.
- Dissolved reagents are stable for at least 2 weeks if stored protected from light and cooled (4°C).
- Dissolved reagents can be stored long-term at -20°C (within expiry). Avoid multiple freeze-thaw cycles.

2. Additional Reagents required

LightCycler[®] Multiplex DNA Master

Cat.-No. 07 339 585 001

LightCycler[®] Multiplex RNA Virus Master

Cat.-No. 06 754 155 001

or Roche LightCycler[®] 480 Probes Master (no instructions included)

Cat.-No. 04 707 494 001

3. Introduction

This product is intended to be used as **Internal Control** to ensure for parameter-negative samples that PCR is working and no inhibitors are present. The amount of target included in the reagent is adjusted to yield a Cp value in the range of 29-31.

The kit 66-0901-96 contains separate target to be used as spiked **Extraction Control** (EC).

4. Description

A 85 bp long fragment from the cloned Phocine herpesvirus (PhHV) sequence target is amplified with specific primers and detected with a LC670 labeled hydrolysis probe.

5. Specification

This assay shall yield a Cp value of about 30 for analytical-PCR negative samples in multiplex PCR.

6. Sample Material and Extraction

Not applicable - depends on the analytical PCR.

7. Material Safety Data (MSDS)

According to OSHA 29CFR1910.1200, Australia [NOHSC:1005, 1008 (1999)] and the EU Directives 67/548/EC and 1999/45/EC any products which do not contain more than 1% of a component classified as hazardous or classified as carcinogenic do not require a Material Safety Data Sheet (MSDS).

Product is not hazardous, not toxic, not IATA-restricted. Product is not from human, animal or plant origin. Product contains synthetic oligonucleotide primers and probes.



8. Instructions for Use

- Instrument programming see document
- Color Compensation see instructions in
- Pipetting instructions multiplex PCR see

ModularDx Programming

40-0320 Universal Color Compensation Hexaplex

ModularDx Multiplex

8.1. Programming Roche 480 Instruments

See the Instrument operator's manual for details. Start programming before preparing the solutions. The protocol consists of three program steps:

- 1: Denaturation: sample denaturation and enzyme activation
- 2: Cycling: PCR-amplification
- 3: Cooling: cooling the instrument

Detection Format 660 Channel

LightCycler® 480 Instrument:

LightCycler® 480 II Instrument:

cobas z 480 Analyzer (open channel):

Set Quant Factor 10, Max Integration time 3 sec

615-670

618-660

610-670

Program Step:	<i>RT Step</i>	<i>Denaturation</i>	<i>Cycling</i>			<i>Cooling</i>
Parameter						
Analysis Mode	<i>None</i>	None	Quantification mode			None
Cycles	1	1	45			1
Target [°C]	55	95	95	60	72	40
Hold [hh:mm:ss]	00:05:00	00:05:00	00:00:05	00:00:15	00:00:15	00:00:30
Ramp Rate [°C/s] 96	4.4	4.4	4.4	2.2	4.4	1.5
Ramp Rate [°C/s] 384	4.6	4.6	4.6	2.4	4.6	2.0
Acquisition Mode	<i>None</i>	None	None	Single	None	None

* 1-Step RT-PCR optional

Table 1

8.2. Experimental Protocol

- **Sample material:** Use aqueous nucleic acid preparations (e.g. 'High Pure PCR Template Preparation Kit').
- **Negative control:** Always run at least one no-template control (NTC) - replace the template NA with water.

For an increased sensitivity use 10 µl sample per 20 µl reaction, in case that inhibition is likely to occur, e.g. extracts obtained from fecal samples, use 5 µl. For 10 µl reactions in 384 well plates use 5 µl /2.5 µl.

8.2.1. Preparation of Parameter-Specific Reagents (PSR, 96 reactions):

One reagent vial with a **green** cap contains all primers and probes to run 96+ LightCycler® reactions.

Add 50 µl PCR-grade water to each reagent vial, mix the solution (vortex) and spin down.

For robotic pipetting the volume can be extended to 55 µl (signals will decrease by 10-20%).

► **Use 0.5 µl** reagent for a 20 µl PCR reaction.

8.2.2. Preparation of the Positive Control

Not applicable - target DNA is included in the reagent.

8.2.3. Preparation of the Reaction Mix

Multiply volumes by the number of reactions plus controls and one reserve and prepare in a cooled tube. The table contains the amounts for a duplex reaction for one analytical parameter (red) and the control:

For use with the Roche LightCycler® Multiplex RNA Virus / DNA Master		
for 5 µl extract	Component	10 µl extract
10.0 µl	Water, PCR-grade (colorless cap, provided with the Roche Master kit)	5.0 µl
0.5 µl	PhHV Control Reaction and	0.5 µl
0.5 µl	PSR Reagent mix additional assay(s) (Multiplex PCR)	0.5 µl
-	For each additional assay reduce the amount of water by 0.5 µl	-
4.0 µl	Roche Master (see Roche manual)	4.0 µl
0.1 µl	Optional RT Enzyme (RNA Master only) - reduce water by 0.1 µl	0.1 µl
15.0 µl	Volume of Reaction Mix	10.0 µl


Table 2

Mix gently, spin down and transfer 15 µl (10 µl) per well.

Add 5 µl (10 µl) of sample or control to each well for a final reaction volume of 20 µl. Seal plate and centrifuge.

Start run

9. Typical Results (Data from LightCycler® 480 II system)



Here:1E3 / reaction

Figure 1

10. Reading the Results

Perform data analysis as described in the operator's manual. For multiplex assays select the color compensation. We recommend using the Second Derivative Maximum method (Automated (F'' max)). View results in the 660 channel. The negative control (NTC) and nergative samples must show a signal.

Analytical PCR (sample)	Channel 660 Control Reaction	Analytical PCR NTC Control	Result
No amplification	Detectable	Negative	Parameter - Negative
Amplification Cp < 37 ⁺	Not relevant	Negative	Parameter - Positive
No amplification	Not detectable	Not relevant	PCR failure Repeat
Amplification signal	Not relevant	Positive	Contamination Repeat

Notes: cobas z 480 Analyzer signal levels are about 50% compared to LightCycler® 480 II results.

+ Recommendation : Define the cut-off 2-4 cycles higher than observed for 10 copies.

11. References

Diagnosing herpesvirus infections by real-time amplification and rapid culture. van Doornum GJ, Guldemeester J, Osterhaus AD, Niesters HG. J Clin Microbiol. 2003 Feb;41(2):576-80

12. Multiplex PCR Compatibility

The PhHV internal control can be combined with up to three (LightCycler® 96), four (cobas z 480) or five (LightCycler® 480 system) analytical assays :

Multiplex PCR and Instrument Compatibility
Color Compensation 40-0320 is mandatory for Multiplex PCR


500	530	580	610	640	660
	Assay 1				PhHV
	Assay 1	Assay 2	Assay 3		PhHV
	Assay 1	Assay 2	Assay 3	Assay 4	PhHV
Assay 5	Assay 1	Assay 2	Assay 3	Assay 4	PhHV

480 II	z 480	LC96	LC2.0	Nano
X	X	X		
X	X	X		
X	X			
X				

Table 3

13. Version History

V140404	Release version	2014-04-14
V140909	Editorial changes	2014-09-09
V150404	LightCycler® Nano removed, measurement at 60°C recommended	2014-04-10
V151001	2015 protocol Multiplex Master, 5/10 µl extract and 60°C acquisition	2015-10-01
V160313	8.2.3 wording	2016-07-07

Certificate of Analysis (CoA)		
Lot n° Expiry :		
Dilution	PSR	passed
Cp range	28-31	
Measured		
Signal level	10-20	
Measured		
Negatives	--	
Note: Cp (crossing point) values collected with pDNA (single target PCR). Fluorescence (FL) levels depend on instrument settings and may vary. The Cp values will vary from instrument to instrument by up to 2 cycles, while the distance between two dilution steps should be relative constant (Δ Cp).		
QC Acceptance Date:		YYYYMMDD
We, the undersigned, certify that the product designated above has been obtained in accordance with the rules of production and quality control.		
Name(s) :		

TIB MOLBIOL Syntheselabor GmbH | Eresburgstr. 22-23 | D-12103 Berlin | Germany
Tel. +49 30 78 79 94 55 | FAX +49 78 79 94 99 | dna@tib-molbiol.de | WWW.TIB-MOLBIOL.COM
Geschäftsführer (CEO): Olfert Landt | Register HRB 93163 B | Registergericht Berlin Charlottenburg

Distributed by Roche - www.lifescience.roche.com